

## Complete Bottle Management

- Fill Scripts Faster
- Automate Perpetual Inventory
- Lower Operating Costs
- Save Space





Drawer Lock provides extra security for controlled meds

MedPickRx lights up to show med location



### Complete Bottle Formulary

Increases...

↑ Fill Accuracy

↑ Inventory Accuracy

↑ Inventory Capacity

And Simultaneously Decreases...

↓ Fill Times


↓ Inventory Stocking Levels

↓ Expired Medications

**MedPickRx** is the latest in modular pharmacy automation. Innovative light directed and barcode verified technology handles the filling process for 100% of your bottle formulary. Perpetual inventory tracking provides powerful reporting capabilities to identify areas for increased savings. MedPickRx can be implemented with little disruption to your pharmacy and results in reduced inventory levels, improved efficiency and increased accuracy.

ORGANIZED BY ACTIVITY • DIRECTED BY LIGHT • VERIFIED BY BARCODE


### BEFORE MedPickRx



## MedPickRx Enhanced Pharmacy

“We were able to add a complete closed-door pharmacy in the same space we previously used for just our retail pharmacy. We now have two pharmacies in the same space!”

### AFTER MedPickRx Installation



- Increased Space**
- 340B
  - Compounds
  - WillCall
  - Workflow
  - Vaccines
  - Consulting
  - Clinic

### BEFORE MedPickRx

- ☀ Constrained growth & profitability
- ☀ Script fill errors (incorrect med or dosage)
- ☀ Time wasted searching for inventory items
- ☀ Only 9% of Oral Solid formulary, 43% of scripts (robotic automation)
- ☀ Expired inventory ranges from 4-8%

### AFTER MedPickRx

- ☀ Reduced labor and space yielding growth potential
- ☀ Technology guided selection and verification for drug and dosage - Results in fewer errors
- ☀ Computerized and optimized inventory reduces fill time by up to 40%
- ☀ Handles 100% of bottle formulary, 92% scripts
- ☀ Nearly eliminates expired medications and on hand inventory reductions

# SOFTWARE

## Scan Order



## 1. Enter Rx

## Interface Order\*



MedPickRx works with or without an interface

User Directed by light to med location



## MedPickRx Technology

- Guides to open bottle first
- Shortest expiration
- Image verification
- Rx detail with status

“We used to have every bottle of a given drug open – all at the same time.”

\* Optional Feature

# SOFTWARE

## 2. Pick Medication, Fill Rx

### Pick Med



- Auto update of perpetual inventory
- Reduced fill time / less walking
- Light direction & barcodes increase accuracy



7 colors for 7 users



“I’ve noticed I can fill a Rx so much faster now. I’ve processed a single Rx in less than 20 seconds. That’s faster than a robot!”

## 3. Return & Stock

### Return & Stock



- Stock randomly - no alphabetizing
- Quick and easy expiration date entry

BEACON						
BEACON EXPIRATION DATE BARCODE SCAN SHEET						
JAN 02	JAN 07	JAN 12	JAN 17	JAN 22	JAN 27	JAN 31
FEB 02	FEB 07	FEB 12	FEB 17	FEB 22	FEB 27	FEB 28
MAR 02	MAR 07	MAR 12	MAR 17	MAR 22	MAR 27	OPEN BOTTLE
APR 02	APR 07	APR 12	APR 17	APR 22	APR 27	EMPTY SLOT
MAY 02	MAY 07	MAY 12	MAY 17	MAY 22	MAY 27	PUR BOTTLE
JUN 02	JUN 07	JUN 12	JUN 17	JUN 22	JUN 27	LOCK BOTTLE

# REPORTS

## MedPickRx Report Generation

Information is valuable - MedPickRx provides powerful reports that help manage your pharmacy's inventory, labor and efficiency.

**Suggested Reorder Report**

Medication	Per #	Max	Suggested Order
CLONAZEPINE	5	10	1
ALPRAZOLAM	1	10	1
AMITRIPTYLINE HCL	1	11	1

**Control Medications Report**

NDC	Medication	Level	Plots	Qty	Barcode
0214809713	ANESTROZOLE 15 mg	Legend	1	20	0214809713
0214809721	ACETAMINOPHEN 325 mg	Legend	1	1	0214809721
0084003704	ALENDRONATE SODIUM 70 mg	Legend	1	75	0084003704

**NDC Usage Report**

NDC	Medication	Barcode	Qty On Hand	Last Use	Usage
0002241530	CLONAZEPINE 5 mg	0002241530	5.00	02/18/2011	90.00
02442300	CLONAZEPINE 20 mg	02442300	30.00	10/12/2010	75.00
0002241530	CLONAZEPINE 5 mg	0002241530	90.00	10/19/2010	200.00

## Innovative Reporting

- Control Medications
- Inventory & Cost
- Medications Dispensed & Remaining
- Product Levels
- Will Expire
- Re-Order
- Returnables
- Unaccessed

“ Since installing Beacon, we have reduced our on-hand inventory from \$237k to \$158k - all in the past six months. ”

Customer Rx Id: **RX253203288**

Patient Name: **9:04 AM\_634351718554375000, AUTOGENRX**

Status	Time	User
Created	3/8/2011 9:04:15 AM	
Dispense_Ready	3/8/2011 9:04:15 AM	admin
Dispense_Active	3/8/2011 9:04:15 AM	admin
Dispense_Ready	3/8/2011 9:07:10 AM	admin
Dispense_Active	3/8/2011 9:33:03 AM	admin
Dispense_Ready	3/8/2011 9:33:32 AM	admin

## Complete Med History

- The Rx
- Who filled it
- Rx fill time
- Filled Quantity
- Inventory Quantity

# STANDARD MODULES

## • 2-BAY HIGH DENSITY MODULE •



9' W x 3.66' D x 6.66' H

The 2-Bay high density unit stores approximately 2,000 - 4,100 items in about 33 square feet.

## • 3-BAY HIGH DENSITY MODULE •



12' W x 3.66' D x 6.66' H

The 3-Bay high density unit stores approximately 3,000 - 6,150 items in less than 40 square feet.

## • 3-BAY WORKSTATION & 1-BAY HIGH DENSITY MODULE •



11' W x 3.66' D x 6.66' H

The 3-Bay Workstation provides a work area for pharmacists and technicians to dispense and check prescriptions with quick access to approximately 2,200 - 3,250 items.

# ADD-ONS

## • CAROUSEL MODULE •



81" Tall  
45" Wide

The Carousel Module stores up to 1,600 items in less than 14 square feet.

- High Density Circular Shape
- Directed by Light Integration
- Unit of Use Storage

## • AVF-70 AUTOMATIC VIAL FILLER •



- Automation for the Top 70 Oral Solids
- Simple, Reliable Design
- Easy Filling and Refilling
- Easy Access Cleaning
- Compatible with All Vials

## • SANYO REFRIGERATION •



- Pharmaceutical Grade
- Precise & Stable Temperature
- Compact
- Reliable

